

Epoch Life Science, Inc.


**P.O. Box 16202
Sugar Land, TX 77496-6202
USA**

Technical Support:

Tel: 832-886-5231
Fax: 832-415-9502
E-mail: support@epochlifescience.com
URL: www.epochlifescience.com



GenCatch™ Plant Genomic DNA Purification Kit



User's Guide for **Plant Genomic DNA Purification**

For Research Use Only

© 2009 by Epoch Life Science, Inc.

www.epochlifescience.com



Overview

GenCatch™ Plant Genomic DNA Purification System provides a simple, fast and cost-effective method to purify genomic DNA from various plant species. It comes with shearing tubes for simple and fast homogenization and filtration of plant tissues. Genomic DNA of predominantly 20-30 kb can be isolated (from 100 mg plant material or 1×10^8 cells) free of contaminants and without phenol/chloroform extraction and ethanol precipitation.

Kit contents:

PX1 Buffer (1), PX2 Buffer (1), PX3 Buffer (1), RNase A(1), WS Buffer (2), Plant Genomic DNA Mini Column (250), Shearing Tube (250), Collection Tube (500) and protocol (1)

Protocol:

<Note>:

- All centrifugation should be done at room temperature.
- Preheat a water bath to 65°C.
- Preheat TE or ddH₂O to 65°C for DNA elution.
- PX1 Buffer and PX3 Buffer may form precipitates, warm at 65°C to dissolve.
- Add 1100 µl of ddH₂O to the RNase A powder tube, vortex to dissolve and store at 4°C.
- Add 180 ml of ethanol (96-100%) to the WS Buffer bottle when first open the bottle.

1. **Grind 100 mg (or less) plant sample under liquid nitrogen to a fine powder and transfer to a new tube.**
Do not allow the sample to thaw, and continue immediately to step 2.
2. **Add 400 µl of PX1 Buffer and 4 µl of RNase A stock solution (100 mg/ml) to the tissue powder and vortex vigorously, then incubate the mixture at 65°C for 10 minutes.**
Do not mix PX1 Buffer and RNase A prior to use. Invert 2-3 times during 65°C incubation.
3. **Add 130 µl of PX2 Buffer to the lysate, vortex, and incubate on ice for 5 minutes.**
4. **Apply lysate to the Shearing tube sitting in a Collection tube and centrifuge at full speed for 2 minutes. Transfer flow-through sample from the Collection tube to a new tube.**
Avoid pipetting any debris or pellet in the collection tube.
5. **Add 0.5 volume of PX3 Buffer and 1 volume of 96-100% ethanol to the clear lysate and mix by pipetting.**
For example: if 450 µl clear lysate collected, add 225 µl PX3 Buffer and 450 µl ethanol.

6. **Apply 650 µl of the ethanol added sample (including any precipitate) from step 5 to a Plant Genomic DNA Mini Column sitting in a Collection tube, close the cap, centrifuge at 8,000 x g (10,000 rpm) for 1 minute, and discard the filtrate.**
If the solution remains above the membrane, centrifuge again at 13,000 rpm.
7. **Repeat step 6 for the rest of the sample.**
8. **Wash the column twice with 0.7 ml of WS Buffer by centrifuging at full speed (13,000 rpm or 10,000 x g) for 30 seconds and discard the filtrate.**
Remember to add 180 ml of ethanol (96-100%) to the WS Buffer bottle when first open the bottle.
9. **Centrifuge at full speed for another 2 minutes to remove traces of WS Buffer.**
10. **Transfer the column to a new 1.5 ml tube, add 200 µl of preheated (65°C) TE or ddH₂O, and centrifuge at full speed for 1 minute to elute DNA.**
11. **Store DNA at -20°C.**